



Product Name :
LecStrip™ M-3 RT Long Wires
Reference :
A1-30R2-2L
Description :

- Set of 3 pre-wired, pre-positioned and pre-gelled ECG electrodes
- Single use/patient
- Repositioning gel
- Radio-transparent

Identification Label :
A1-30R2-2L

**ADULT
M-3
RADIO TRANSPARENT**
Limb : I, II, and III

Use :

- Adult ECG monitoring
- Up to 3 leads
- Up to 72 hours
- X-ray

Compatibility :

- LecStrip specific ECG trunk cables:
ref. LMXX-3 (3 leads monitoring)
"XX" indicates variant number ECG brand compatibility and are available on catalogue

Packaging :

- Type
 - Individual aluminium foil pouch per set (3 electrodes)
 - Box of 40 pouches
 - Carton of 12 boxes
- Information displayed on box and pouch
 - Traceability label (lot No, Manufacturing date, Expiry date using UDI number GS1 format)
 - Identification label (see opposite)
 - CE mark (see below)

Sensor Specification :

- Hydrogel 33x33mm, thickness 0.85mm +/- 0.05mm
- Non-woven fabric back support
- Adhesivity 75 gr/cm minimum
- Impedance 200 Ohms max. à 10 Hz
- Volume resistivity 1000 Ohm-cm maximum
- Bio-compatibility (tests cytotoxicity / irritation / hypersensitivity)

Wire Specification :

- Medical grade UL-PVC insulation
- Carbon conductor (radio-transparent)
- Impedance 145 Ohms/m +/- 10 Ohm
- LecStrip Connector (molex type) with protective comfort cap
- Withstand defibrillation shocks

CE Mark
Class I
 Annex IX directive 93/42 modified by
 Directive 2007/47/CE

Registering Date
 14/11/2016

Compliant With Standards :

- IEC 60601-1
- ANSI/AAMI EC-12:2000/(R)2005
- ANSI/AAMI EC-53:1995/(R)2001
- IEC 60601-2-25:1993+A1:1999, 51.101.1 (report WT093051)

Manufacturer :

 ASFALIS S.A.S.
 14 résidence Ker An Kloareg
 29170 Fouesnant
 FRANCE
 Email: contact@asfal.fr