



TECHNICAL DATA SHEET FOR 4-PLY H.C.L FACE MASK

Art No	Item	Style	Standard	Country of Origin
	H.C.L face mask	4-ply ear-loop	Type IIR	Viet Nam
		Description: Rectangular face masks with a shapeable nosepiece and two earloops present, one on each side, in order to hold mask in place		
		Mask body	Layer 1	Spunbond Polypropylene
			Layer 2	Spunbond Polypropylene
			Layer 3	Meltblown
			Layer 4	Spunbond Polypropylene
Nose-piece	Aluminum / Plastic covered iron			
Earband	Elastic earloop creates the comfortable feeling			

Dimension	Adult		Kid
Body size	Length	175 mm	150 mm
	Width	95 mm	95 mm
Length of Earloop	180 mm	145 mm	
Length of Nose-piece	115 mm	90 mm	

	Color	Blue, White, Black		
	Standard	Type IIR	BFE(Bacterial filtration efficiency) $\geq 99\%$	
			Delta P(Differential pressure)(Pa/cm ²) < 29.4	
			Delta P(Differential pressure)(mm H ₂ O/cm ²) < 4	
			SBP (synthetic blood penetration) 120mmHg ≥ 16	
Packing	Package	50pieces/box, 50boxes/carton		
	Box size	18.5*10*12 cm		
	Carton size	52*40*64 cm		

Usage instruction	Open box/bag and check face mask for damage and signs of pollution. When wearing the face mask, keep the nose-clip side up on your nose, press and adjust the nose-clip to make sure it fits the nose area.			
	 FRONT	 BACK	 OPEN THE LEFT EARLOOP ENSURE YOU OPEN WITH CLEAN HANDS	 OPEN BOTH AND WEAR IT

Usage suggestion	Continuously wear the face mask - do not use for more than four hours. Single use only - do not re-use
Storage suggestion	Store in a clean, dry environment and avoid excessive heat

Contact Us

Head Office: Hieu Linh DTTM Joint Stock Company
Address: No 2, Lane 94/17 Ngoc Khanh Str, Giang Vo Ward, Ba Dinh District, Ha Noi, Viet Nam
EU Office: Vfood Aps
Address : Lokesalle 69 8700 Horsens, DK.
Tel: +45 42 668 668
Email: vd@vfood.dk