

Data Sheet

1/ Product Name: Disposable Injection Needle

2/ Reference Number: VDK-IN-23-230-2304-A

VDK-IN-23-230-2305-A

VDK-IN-23-230-2303-A

3/ Intended Use: It is used with endoscope for the injection treatment of esophageal and gastric varices, and the injection of submucosa in GI tract.

4/ Features: It has high-flow lumen for smooth injection with light force.

Shielded needles designed to prevent damage to working channel.

5/CE classification: Class IIa

6/ Materials (Table1)

Table1

Component	Material	Standard
Needle	X5CrNi18-10	EN 10088-1-2005
Inner Tube	PTFE	ISO 13000-2005
Injection Handle	ABS	ISO 2580-1-2002
Boosting Tube	X5CrNi18-10	EN 10088-1-2005
Front Handle	ABS	ISO 2580-1-2002
Protective Sleeve	PEBAX	ISO 1874-2010
Outer Tube	PTFE	ISO 13000-2005
Guiding Head	X8CrNiS18-9	EN 10088-1-2005

