

**Product Name :**

**ECG Interface Cable for LecStrip™ Electrodes Lead sets  
5 Leads**

**Reference :**


**LM050-5F**

**Description of the device :**

Dimensions : 3,4 m  
 ECG connector to yoke: 0,7 m  
 Yoke to LecStrip connector: 2,7 m

Weight: 210 gr  
 Cable material: TPU  
 Yoke / Connectors material: Medical Grade PVC  
 Presence of latex: No  
 Presence of DHP: No

For further details on use and maintenance see Instructions For Use


**Compatibility with :**

Electrodes: 5 leadwires LecStrip electrodes; ref. N1-50R2-1 / N1-50R2-1L / A1-50R2-1 / A1-50R1-1L

ECG Machine: Philips series MP, MX, VM, C, Merlin, MMS

**Packaging :**

1 cable per bag  
 Unit of order, Qty: 1 bag  
 Minimum quantity of delivery: 1

**Regulatory :**

CE Mark: Yes

Applicable EU regulation: MDR 2017/145

Medical Device Class: I

Classification rules: Annex VIII

Number of certifying body: N/A for Class I

Traceability: Laser marked UDI number with serial number on cable


**Standards :**

IEC 60601-1, IEC 60601-2-25 (report WT093051)  
 ANSI/AAMI EC-12, EC-53  
 Bio-compatibility (ISO tests cytotoxicity / irritation / hypersensitivity)

**Guarantee :**

The product is guaranteed for 1 year against manufacturing defects.  
 The defects and the deteriorations caused by the natural wear or by an outside accident, or by the modification of the product not planned or not specified by the seller are excluded from the guarantee.

**Manufacturer :**

ASFALIS S.A.S.  
 14 résidence Ker An Kloareg, 29170 Fouesnant, FRANCE  
 Email: contact@asfalis.fr