

Product Name :

**ECG Interface Cable for LecStrip™ Electrodes Lead sets
3 Leads**


Reference :

LM050-3F

Description of the device :

Dimensions :	3,4 m
	ECG connector to yoke: 0,7 m
	Yoke to LecStrip connector: 2,7 m
Weight:	210 gr
Cable material:	TPU
Yoke / Connectors material:	Medical Grade PVC
Presence of latex:	No
Presence of DHP:	No

For further details on use and maintenance see Instructions For Use


Compatibility with :

Electrodes:	3 leadwires LecStrip electrodes; ref. N1-30R2-1 / N1-30R2-1L / A1-30R2-1 / A1-30R1-1L / A1-40R2-1
ECG Machine:	Philips series MP, MX, VM, C, Merlin, MMS

Packaging :

1 cable per bag	
Unit of order, Qty:	1 bag
Minimum quantity of delivery:	1

Regulatory :

CE Mark:	Yes
Applicable EU regulation:	MDR 2017/145
Medical Device Class:	I
Classification rules:	Annex VIII
Number of certifying body:	N/A for Class I
Traceability:	Laser marked UDI number with serial number on cable

Standards :

IEC 60601-1, IEC 60601-2-25 (report WT093051)
 ANSI/AAMI EC-12, EC-53
 Bio-compatibility (ISO tests cytotoxicity / irritation / hypersensitivity)

Guarantee :

The product is guaranteed for 1 year against manufacturing defects.
 The defects and the deteriorations caused by the natural wear or by an outside accident, or by the modification of the product not planned or not specified by the seller, are excluded from the guarantee.

Manufacturer :

ASFALIS S.A.S.
 14 résidence Ker An Kloareg, 29170 Fouesnant, FRANCE
 Email: contact@asfalis.fr